



VISTA  
LABORATORY  
SERVICES

### Recommended Antibody Pairing - 2021

| <u>Target Antigen</u> | <u>Capture Antibody</u> | <u>Detector Antibody</u> |
|-----------------------|-------------------------|--------------------------|
| Falciparum pLDH       | 17e4                    | 12g1 or 6F1              |
| Vivax pLDH            | 3h8                     | 12g1 or 6F1              |
|                       | 13h1                    | 12g1 or 6F1              |
|                       | 1102                    | 1201                     |
| Ovale pLDH            | 10d1                    | 12g1 or 6F1              |
| Pan Malaria pLDH      | 19g7 +10d1              | 12g1 or 6F1              |
|                       | IB47                    | DH60                     |
| Falciparum HRPII      | 6F12                    | 9F2                      |
|                       | 4F4                     | 1F12                     |
|                       | 4H44                    | 4H45                     |
| HBsAg                 | pAb-HBs                 | mAb-HBs                  |