



VISTA  
LABORATORY  
SERVICES

## Recommended Antibody Pairing - 2022

<u>Target Antigen</u>	<u>Capture Antibody</u>	<u>Detector Antibody</u>
Falciparum pLDH	17e4	12g1 or 6F1
Vivax pLDH	3h8	12g1 or 6F1
	13h1	12g1 or 6F1
	1102	1201
Ovale pLDH	10d1	12g1 or 6F1
Pan Malaria pLDH	19g7 +10d1	12g1 or 6F1
	IB47	DH60
Falciparum HRPII	6F12	9F2
	4F4	1F12
	4H44	4H45
HBsAg	pAb-HBs	mAb-HBs