



**VISTA
LABORATORY
SERVICES**

Recommended Antibody Pairing - 2022

<u>Target Antigen</u>	<u>Capture Antibody</u>	<u>Detector Antibody</u>
Falciparum pLDH	17e4	4c4, 6F1 or 12g1
Vivax pLDH	3h8	4c4, 6F1 or 12g1
	13h1	4c4, 6F1 or 12g1
	1102	1201
Ovale pLDH	10d1	4c4, 6F1 or 12g1
Pan Malaria pLDH	19g7 +10d1	4c4, 6F1 or 12g1
	IB47	4c4 or DH60
Falciparum HRPII	6F12	9F2
	4F4	1F12
	4H44	4H45
HBsAg	pAb-HBs	mAb-HBs