



VISTA  
LABORATORY  
SERVICES

### Recommended Antibody Pairing - 2023

<u>Target Antigen</u>	<u>Capture Antibody</u>	<u>Detector Antibody</u>
Falciparum pLDH	4C8 or 17e4	4c4, 6F1 or DH60
Vivax pLDH	3h8	4c4, DH60, 3h8
	13h1	4c4, 6F1 or 12g1
	1102	1201
Ovale pLDH	10d1	4c4, 6F1 or 12g1
Pan Malaria pLDH	19g7 +10d1	4c4 or DH60
	IB47	DH60
Falciparum HRPII	4H48	4H46
	4F4	1F12
	6F12	9F2
HBsAg	pAb-HBs	mAb-HBs