



VISTA
LABORATORY
SERVICES

Recommended Antibody Pairing - 2024

<u>Target Antigen</u>	<u>Capture Antibody</u>	<u>Detector Antibody</u>
Falciparum pLDH	4C8* or 17e4	4c4*, 6F1 or DH60
Vivax pLDH	3h8*	4c4, DH60, 3h8*
	13h1	4c4*, 6F1 or 12g1
	1102	1201
Ovale pLDH	10d1	4c4, 6F1 or 12g1*
Pan Malaria pLDH	19g7*+10d1	4c4 or DH60*
	IB47*	DH60*
Falciparum HRPII	4H48*	4H46* - 4H45*
	4F4	1F12
	6F12	9F2
HBsAg	pAb-HBs	mAb-HBs

* = Our Recommended Pairing